Intergrip

en

Dear Customer, We appreciate the confidence placed on a product EMO. Choosing a brace Intergrip brand you have purchased a quality product and high level medical doctor. Read these instructions carefully. If you have any questions please contact your doctor or your specialist.

REGULATIONS

The medical device complies with Regulation (EU) 2017 /745 of the European Parliament and of the Council on medical devices and the respective national regulations. Class I medical device, non-sterile, without measurement function.

DESCRIPTION AND FEATURES

The IF702 Intergrip is an immobilizing orthosis of the interphalangeal joints. Soft, light and padded fabric. Incorporates conformable aluminum splint. It incorporates removable wrist pad. Velcro fasteners in different positions. Velcro adjustable.

INDICATIONS

Joint trauma Dislocations and sprains. Finger fractures. Capsulitis. Bursitis. Dupuytren contracture. Post-surgical treatment.

EFFECTS

Immobilization of interphalangeal and metacarpophalangeal joints. Protection after surgery. Pain relief.

INSTRUCTIONS FOR USE AND PLACEMENT

The product must be adapted by the Orthopedic Technician or by the physician individually for each patient. The correct choice of size is essential to achieve maximum treatment efficiency, as well as prolong the life of the garment. Check the conformation of the strap and adapt to morphology if necessary before placing it on the patient. Open all velcro closures and place the finger or fingers to be treated on the orthosis. Next, fix the finger, hand and wrist tapes in this order.

Once immobilized, check that the straps do not press excessively and that the treatment is comfortable for the patient. The patient can be removed and put on the brace without hardly mobilizing the limb, facilitating rehabilitation exercises and personal hygiene if these actions have been indicated by the technician. The orthosis is served with a pad to place, if necessary, on the groove of the ventral wrist face.

PRECAUTIONS

Follow the manufacturer's directions for proper orthosis placement. Follow the instructions or recommendations of the prescriber or healthcare personnel. For a personalized adaptation, consult a specialist. If there is ongoing pain, it is advisable to contact a healthcare professional. The used material is hypoallergenic but if there is any allergic reaction or irritation remove the product and consult with the prescribing doctor. The user should position itself in such a way that the compression is not excessive. In case of doubt about the size, choose the bigger size. Only use

during active periods.

Remove during the night. The correct choice of the appropriate size for each patient is essential to obtain the highest degree of therapeutic efficacy. The adjustments of the orthotic should be made individually. The orthotic should be used for the purpose it is prescribed for. Do not expose the orthotic to a flame, heat source or high temperatures.

CLEANING AND MAINTENANCE

Hand wash with warm water (max. 30°) and neutral soap. During washing, if the garment has Velcro, they must stick together or keep them closed. Dry at room temperature. Do not iron the garment shop. Contact with ointments, creams, alcohol and solvents can damage the garment.

QUALITY CERTIFICATES

All the materials used have been tested, comply with the applicable European regulations, are of high quality and offer maximum safety and comfort during use. Öko-Tex Standard 100 Certificate: Ecological textile label granted after the creation of a catalog of analysis criteria that guarantees the absence of substances harmful to health. It certifies the absence of harmful substances in all the materials with which a garment is made. substances in all the materials with which a garment is made.

WARRANTY

Especialidades Médico Ortopédicas s.l., guarantees that all its products, as long as they have not been manipulated or altered, always are within the framework of the legal terms. Products whose characteristics have been altered due to misuse are not included in this warranty.

NOTICE

To dispose of the product, use the solid waste container. For the disposal of the container, depending on its material, use the paper / cardboard container or the plastic container. In any case, you must strictly comply with the regulations of the country. Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and / or the patient are established.

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