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ORTHOS

inclusive solutions
WALKING AIDS

Instructions Manual WALKER WITH WHEELS

ORTHOS XXI products are equipped with studied equipments in order to serve de patients, handicaps and elders, with multiple handicaps; its utilization must be prescribed by health professionals or authorized salesmen who are familiar with the products, what they are for and the instructions manual.

Our wide experience and manufacture technology allow us to adjust the products according to the physical handicaps of the final user. The walkers with wheels are Class I medical devices according to the communitary directive 2007/47/EC of September 5th and tested according to the ISO 11199-2:2005.

Thank you for choosing ORTHOS XXI products.

MAIN FEATURES

The walkers with wheels are external use walking aids which replace partly or fully a lower limb and which are to be used by people temporarily or permanently handicapped, patients or elders in their go about on foot.

The walkers with wheels are made of an aluminium tube with a handgrip adjustable in height by metallic pins, with rubber tips in the end touching the walking ground. The device has a finishing in eletrostatic painting in poliester powder that possesses a good mechanic resistance to impact and to UV rays.

The walkers with wheels are made of an aluminium tube with adjustable in height handgrips by means of strings with 2 pins, with 2 rubber tips and 2 wheels at the end touching the walking ground.

TO ASSURE THE LONGEVITY AND EFFICIENCY OF YOUR WALKER WITH WHEELS YOU MUST:

- Be aware that the walker must not be inappropriately used (hitting doors, throwing objects, impacts or strong hits).
- Replace the tips every time its appearance shows wearing off or deterioration.

WAY OF USE

- Adjust the height of your walker correctly accordingly to the user's height.
- Keep in mind that when the aid is used away from the body or badly supported it may cause accidents due to unbalance or falling.
- The walking aid must be held by the handle grips placed on the sides.
- The walking should be in small steps and using the body movement.

SPECIAL CARES

- Be aware that wet flat surfaces or slippery grounds which may cause the instability of the walking and the user's unbalance and/or fall.
- The soft surfaces are highly dangerous because the tip or the wheels may plunge and cause unbalance to the user.
- The irregular grounds (stones, grids on the sidewalks, small steps) may strain the tips or obstruct the walker which can cause unbalance

and/or falls. Therefore walking stick should be well raised and supported when going through stony and bumpy grounds.

- Pay special attention when walking around through grounds in contact with chemical products because it can corrode prematurely the tips and wheels.

PRESCRIPTION SKILLS

Only the medical staff such as doctors, physiatrists, orthopaedists, nurses or other health professionals have the competence to verify the patients, the physically disable individuals and the elders ability to use these walking aids. In alternative and when the medical staff is not available, the pharmacy auxiliary staff, the official orthopaedic materials' dealer, the assistants of home for elders or the relatives of the physically disabled, after reading the instructions, may advise about the adequate height adjustment and check the condition of the tips and wheels and replacing them when they are wore off or damaged.

The use of walkers with wheels should occur by necessity, not by the device itself but by the disabilities of its user. They are a walking assistance device. Those who use them are less dependent, freer, and able to do their daily walks, go to work, leave their houses, have fun and live in society.

Model	Maximum dimensions of the base (cm)	Height minimum/maximum (cm)	Maximum height of user (cm)	Maximum weight of user (kg)
Adult	64x49	82 / 92	180	100
Child	55x40	74 / 85	150	100

EN ISO 11199-2:2005

NOTE: There is also an bariatric walker that can support a maximum user weight of 300kg.

MAINTENANCE OF THE DEVICE

The device must be clean when necessary with a low abrasive cloth soaked in a neutral detergent diluted in water. Dry the device with a soft cloth. The device can be disinfected with alcohol and derivatives, since they are diluted in water. The thermoplastic components should be inspected regularly and replaced when they show wear or fracture.

Maintenance and repair made by manufacturer

If the manufacturer's intervention is required for maintenance or repair, contact the dealer.

MANUFACTURER WARRANTY

This product has two years warranty after the purchase's invoice date, covering production imperfections. The warranty does not cover the wearing out produced by a normal utilization of the device, as well as defects and damages related with improper manipulation or use or inadequate maintenance.

The transport costs for the repair are of the clients' responsibility.

March 2018

Orthos XXI Unipessoal, Lda.
Rua de Santa Leocádia Nº 2735
4809-012 Stª Leocádia de Briteiros
Guimarães - Braga - PORTUGAL
Tel.: (+351) 253470270 - Fax.: (+351) 253470278
www.orthosxxi.com info@orthosxxi.com